

Microtomy oil

Catalog number	Unit size
08-1720	2 x 100 ml
08-1721	1 x 100 ml

Packaging Primary container: neutral bottle in PE with resealable spout designed to an easy dispensing. Useful capacity 100 ml.
Secondary container: carton box.

Wear, water, alcohol and solvents resistant PVC label. Scratchproof ink resistant to water and alcohol.

Expected aim Product for the preparation of cyto-histological samples for optical microscopy.

Specifications

Density at 20°C	kg/m ³	862
Viscosity at 40°C	cSt	37
Flammability	°C	206
Freezing	°C	-15
Appearance		clear

Application Oil for microtomy

Components

Components	CAS	CE	Index
-	-	-	-

Warning and precaution The product must be used exclusively by specialized technical operators.
Carefully read the information on the classification of dangerous substances on the label. Always refer to the safety data sheet where are available the information on the risks presented by the mixture, the precautionary measures during use, the measures first aid and the intervention in the event of accidental release.

Do not use if the primary container is damaged.

Storage Store the preparation at room temperature. Keep the containers tightly closed.

Stability After the first opening, the product is usable until the expiry date, if correctly stored. Product validity: 5 years.

Disposal Hazardous preparation: observe all state and local environmental regulations regarding waste disposal.

Date of issue: May 2013