

Protect every life story™ in your lab.

TMD™ Tissue Marking Dyes



✓ Expanded Selection of Packaging Sizes and Configurations

We've expanded the packaging options to give you more choices in bottle sizes - so you can use exactly what you need and what works best for your lab. Easy-to-use sizes include 20 ml., 2 oz., and 8 oz. refills. Colors are available individually and in 7-color kits.

✓ New High-Precision Dispenser Cap

Gain new levels of accuracy and control when marking tissue. TMD's precision dispenser cap enables PA's to dispense a precise amount of dye in the exact area they want it, minimizing overflow and wasted dye. The low viscosity allows for a thin, even coating of pigment, ensuring superior adherence and permanent marking of tissue throughout the entire process.

✓ No Mordant or Fixative Required

TMD enables you to apply and dry in one easy step. Our unique formulation includes fixative with the dye in the same bottle. With TMD, there is no need for a separate fixative to be applied after the dye - which is one less step for the PA's.

✓ Quick Drying

TMD requires only 30 - 60 seconds to dry. This advanced formula applies directly, dries quickly, and marks permanently.

✓ Consistent Quality

With TMD, you will be assured of consistent quality throughout the entire package, from when the bottle is first opened to when the last of the dye is dispensed.